

ELISA 2.0 Inhalable Allergens - Validated Performance Parameters

Allergen:	HOUSE DUST MITE							STORAGE	ANIMAL						
	Der p 1	Der p 1a	Der p 23	Der p 2	Der f 1	Der f 2	Blo t 5	MITES	Tyr p 2	Fel d 1	Fel d 4	Can f 1	Mus m 1	Ory c 3	Rat n 1
Linearity (R²)¹	1.0	1.0	0.998	0.999	0.999	0.999	0.999	0.999	0.999	0.999	1.0	1.0	1.0	0.999	1.0
Range (ng/ml)²	100-0.78	50-1.56	50-0.10	12.5-0.10	25-0.39	62.5-0.98	7.5-0.12	250-7.81	25-0.20	5-0.04	25-0.39	25-0.10	20-0.16	50-0.39	
Limit of Quantification³															
<i>LLOQ (ng/ml)^{3a}</i>	0.39-1.56	0.39-1.56	0.391-56	0.10-0.40	0.20-0.39	0.49-1.95	0.12-0.23	1.95-7.80	0.20-0.39	0.02-0.31	0.39	0.10-0.20	0.08-0.31	0.20-0.39	
<i>ULOQ (ng/ml)^{3b}</i>	50-25	50-6.25	25-12.5	25-6.25	25-12.5	250-62.5	30-3.75	500-250	25-12.5	5	25	25-12.5	20-5	50-25	
Accuracy (% Recovery)⁴															
<i>Intra-assay (n=9)^{4a}</i>	98-117%	82-129%	87-129%	83-132%	88-118%	73-118%	84-130%	84-132%	96-129%	93-106%	83-103%	91-113%	92-126%	85-129%	
<i>Inter-assay (n=54)^{4b}</i>	103%	107%	101%	104%	106%	88%	105%	108%	113%	97%	91%	103%	104%	102%	
Precision (%CV)⁵															
<i>Intra-assay (n=9)^{5a}</i>	4-13%	1-23%	2-25%	4-14%	5-15%	1-12%	1-16%	3-14%	4-13%	4-8%	5-18%	5-9%	2-14%	6-12%	
<i>Inter-assay (n=54)^{5b}</i>	8%	11%	8%	7%	10%	7%	8%	6%	9%	6%	9%	8%	6%	10%	

Allergen:	COCKROACH					POLLENS					MOLDS		CANNABIS
	Bla g 1	Bla g 2	Bla g 5	Per a 7	Bet v 1 EP	Phl p 5	Amb a 1	Lol p 1	Cry j 1	Alt a 1	Asp f 1	Can s 3	
Linearity (R²)¹	1.0	0.999	0.999	0.999	0.999	1.0	1.0	0.994	1.0	1.0	1.0	1.0	
Range (ng/ml)²	50-0.39	100-0.39	125-1.95	12.5-0.20	50-0.39	250-0.98	100-0.78	125-7.81	50-0.78	25-0.10	40-0.31	62.5-0.49	
Limit of Quantification³													
<i>LLOQ (ng/ml)^{3a}</i>	0.20-0.78	0.39-1.56	0.98-7.81	0.19-0.39	0.39	1.95-0.98	1.56-0.78	15.63-3.91	0.78-0.39	0.40-0.10	0.16-0.63	0.24-0.49	
<i>ULOQ (ng/ml)^{3b}</i>	50-25	50-25	250-125	50-6.25	100-50	250-62.5	200-25	250-125	50-25	25-6.25	40-20	62.5	
Accuracy (% Recovery)⁴													
<i>Intra-assay (n=9)^{4a}</i>	71-108%	92-113%	77-129%	70-121%	82-120%	85-120%	90-108%	81-112%	94-106%	90-115%	78-115%	88-127%	
<i>Inter-assay (n=54)^{4b}</i>	94%	101%	102%	99%	100%	103%	100%	96%	100%	99%	95%	106%	
Precision (%CV)⁵													
<i>Intra-assay (n=9)^{5a}</i>	3-9%	5-15%	3-16%	2-15%	1-9%	7-10%	6-12%	7-16%	6-16%	3-11%	4-22%	2-19%	
<i>Inter-assay (n=54)^{5b}</i>	6%	10%	8%	6%	7%	9%	8%	12%	10%	6%	12%	11%	

ELISA 2.0 - Validated Performance Parameters

Food Allergens:

Allergen:	Ana o 3	Ara h 1	Ara h 2	Ara h 3	Ara h 6	Ara h 8	Api g 1	Bos d 5 (Native)	Bos d 11	Cor a 9
Linearity (R²)¹	1.0	1.0	0.999	0.999	1.0	1.0	0.999	1.0	1.0	0.998
Range (ng/ml)²	20-0.16	1000-31.25	125-0.98	62.5-0.49	25.0-0.05	12.5-0.20	100-1.56	12.5-0.010	500-7.81	50-0.39
Limit of Quantification³										
<i>LLOQ (ng/ml)^{3a}</i>	0.16	15.63-31.25	0.49-3.91	0.49-1.95	0.05-0.20	0.20-1.56	1.56-6.25	0.10-0.78	1.95-7.81	0.39-0.78
<i>ULOQ (ng/ml)^{3b}</i>	20-5	1000	250-31.25	62.5-31.25	25-12.5	25-6.25	200-50	25-12.5	500-125	100-12.5
Accuracy (% Recovery)⁴										
<i>Intra-assay (n=9)^{4a}</i>	85-110%	77-128%	80-116%	74-127%	105-113%	100-109%	81-130%	84-120%	75-128%	75-112%
<i>Inter-assay (n=36)^{4b}</i>	99%	109%	101%	97%	108%	102%	108%	103%	94%	93%
Precision (%CV)⁵										
<i>Intra-assay (n=9)^{5a}</i>	4-10%	2-19%	7-14%	3-11%	7-13%	3-15%	3-11%	4-9%	2-16%	5-15%
<i>Inter-assay (n=36)^{5b}</i>	8%	8%	10%	5%	8%	8%	7%	6%	7%	10%

Allergen:	Cyp c 1	Gal d 1	Gal d 2	Gly m 5	Jug r 1	Pru du 6	Ses i 1	Shrimp Tropomyosin	Sin a 1
Linearity (R²)¹	0.997	1.0	1.0	1.0	0.997	0.999	0.999	0.999	0.999
Range (ng/ml)²	62.5-0.98	250-7.81	50-0.39	250-1.95	25-0.39	125-0.49	25-0.20	25-0.20	15-0.06
Limit of Quantification³									
<i>LLOQ (ng/ml)^{3a}</i>	0.98-3.90	3.90-7.80	0.39-0.78	1.95-3.91	0.40-0.78	0.49-3.91	0.10-0.40	0.10-0.39	0.03-0.06
<i>ULOQ (ng/ml)^{3b}</i>	125-31.25	250	100-25	250-62.5	50-6.25	125-31.25	25-12.5	25-6.25	15-7.5
Accuracy (% Recovery)⁴									
<i>Intra-assay (n=9)^{4a}</i>	71-120%	89-101%	91-132%	93-119%	86-119%	92-115%	71-131%	95-115%	70-132%
<i>Inter-assay (n=36)^{4b}</i>	86%	95%	104%	105%	108%	102%	109%	108%	96%
Precision (%CV)⁵									
<i>Intra-assay (n=9)^{5a}</i>	2-16%	8-13%	5-16%	5-11%	6-12%	4-15%	3-27%	4-14%	3-16%
<i>Inter-assay (n=36)^{5b}</i>	7%	9%	11%	8%	9%	7%	10%	7%	8%

Validation Notes:

1. **Linearity** is the mean R² of six ELISA plates for control curves generated using 4-parameter logistic fit.
2. **Range** is the average usable range of control curves from six ELISA plates that have a value of 70-130% of the expected concentration, with %CV < 15 between duplicate points.
3. **Limit of Quantification**
 - 3a. LLOQ - The lowest concentration points of six control curves with a recovery of 70-130% and %CV <15, expressed as a range.
 - 3b. ULOQ - The highest concentration points of six control curves with a recovery of 70-130% and %CV <15, expressed as a range.
4. **Accuracy**
 - 4a. Intra-assay - The range of average percent recovery of samples A, B, and C run in triplicate from six ELISA plates (n=9).
 - 4b. Inter-assay - The overall average percent recovery of samples A, B, and C run in triplicate from six ELISA plates (n=54).
5. **Precision**
 - 5a. Intra-assay - The range of average percent coefficient of variation of samples A, B, and C run in triplicate from six ELISA plates (n=9).
 - 5b. Inter-assay - The overall average percent coefficient of variation of samples A, B, and C run in triplicate from six ELISA plates (n=54).