

Product Code: LTN-GD1-1

Gel Image:

LoTox™ nGal d 1

Allergen: nGal d 1 (Trypsin Inhibitor)

Lot No: XXXXX

Source: Partially purified ovomucoid

Mol. Wt: 28 kD

Purification: Ion exchange and size exclusion chromatography.
Purity >95% by SDS-PAGE.

Concentration: See product insert.

Formulation: Preservative and carrier-free in endotoxin free phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert.

Notes: (1) Ovomucoid runs anomalously on SDS-PAGE at 30-45 kDa (3).
(2) A LoTox™ product, Endotoxin < 0.03 EU/µg.



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com