

Product Code: RP-AO3-1

Allergen: rAna o 3 (*Anacardium occidentale* allergen 3.0101)

Lot No: XXXXX

Source: *Pichia pastoris*

Mol. Wt: 15 kD

Purification: Purified from *Pichia pastoris* culture by His-Trap chromatography. Purity > 95 % by SDS-PAGE.

Concentration: 1.0 mg/ml by OD 280. Extinction coefficient = 0.49.

Formulation: Preservative and carrier-free in 20 mM Tris, 0.15 M NaCl, pH 8.2. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: XXXXX

Notes: (1) rAna o 3 has a C-terminal 6xHis-tag

Gel Image:

rAna o 3



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com