

Product Code: RP-CA14-1

Allergen: rCor a 14 (*Corylus avellana* allergen 14.0101)

Lot No: XXXXX

Source: *Pichia pastoris*

Mol. Wt: 15 kD

Purification: Purified from *Pichia pastoris* culture by multi-step chromatography. Purity > 95 % by SDS-PAGE, run under reducing conditions.

Concentration: 1.4 mg/ml by OD 280 nm; Extinction coefficient = 0.39.


Formulation: Preservative and carrier-free in 1X PBS, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: XXXXX

Notes: (1) rCor a 14 has a C-terminal 6xHis-tag.

Gel Image:
rCor a 14



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com