

Product Code: NA-BD4-1

Allergen: nBos d 4 (a-lactalbumin)

Lot No: XXXX

Source: Partially purified bovine milk

Mol. Wt: 14 kD

Purification: HPLC. Purity > 95 % by SDS-PAGE.

Concentration: See Product Insert

Formulation: Preservative and carrier-free in endotoxin free phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See Product Insert

Gel Image:

nBos d 4



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com