

Product Code: NA-CA9-1

Allergen: nCor a 9 (*Corylus avellana* allergen 9)

Lot No: XXXXX

Source: 11S globulin from hazelnut.

Mol. Wt: Multiple subunits ranging from 30 to 59 kD

Purification: Purified from hazelnut extract by multi-step HPLC.
Purity > 95 % by SDS-PAGE.

Concentration: 1.2 mg/ml by OD 280 nm; Extinction coefficient = 0.826

Formulation: Preservative and carrier-free in 1M NaCl, PBS pH 7.40.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: XXXXX

Notes: (1) 11S globulins are synthesized as single peptide chain precursors that are cleaved post-translationally.
(2) Solubility of Cor a 9 is dependent on the presence of at least 1M NaCl.

Gel Image:

nCor a 9



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com