

Product Code: RE-DP10-1

Allergen: rDer p 10 (*Dermatophagoides pteronyssinus*
allergen 10.0101)

Lot No: XXXXX

Source: *E. coli*

Mol. Wt: 34 kD

Purification: Multi-step chromatographic purification. Purity > 95% by SDS-PAGE.

Concentration: 1.0 mg/ml based on BCA Protein Assay.
(Thermo Scientific, Rockford, IL)

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: XXXXX

Notes: (1) rDer p 10 has a C-terminal 8xHis-tag.
(2) Under non-reducing conditions a ~68 kD Der p 10 dimer is visible on SDS-PAGE.

Gel Image:

rDer p 10



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com