

Product Code: NA-LP1-1

Allergen: nLol p 1 (*Lolium perenne* allergen 1)

Lot No: XXXX

Source: Rye grass pollen (*Lolium perenne*)

Mol. Wt: 27 kD

Purification: From Rye grass pollen extract by multi-step chromatography.
Purity > 90 % by SDS-PAGE

Concentration: See product insert.

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Notes: (1) Under non-reducing conditions, a dimer is present at ~55kD in addition to the Lol p 1 monomer at ~27kD.

Gel Image:
nLol p 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com