

**Product Code:** LTN-GD5-1

**Gel Image:**

LoTox™ nGal d 5

**Allergen:** nGal d 5 (*Gallus domesticus* allergen 5)

**Lot No:** XXXXX

**Source:**

**Mol. Wt:** Chicken Serum Albumin  
67 kD



**Purification:** From Chicken Serum Albumin extract by affinity chromatography.  
Purity >95% by SDS-PAGE.

**Concentration:** 1.3 mg/ml by BCA Protein Assay.

**Formulation:** Preservative and carrier-free in endotoxin free phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** XXXXX; When stored at -20°C.

**Notes:** (1) A LoTox™ product, Endotoxin < 0.03 EU/µg.

For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

[www.inbio.com](http://www.inbio.com)

700 Harris Street  
Charlottesville, VA  
22903  
United States

(434) 984-2304  
[mail@inbio.com](mailto:mail@inbio.com)

Vision Court  
Caxton Place  
Cardiff, Wales CF23 8HA  
United Kingdom

+44 (0) 2921 674 640  
[Info@inbio.com](mailto:Info@inbio.com)