

**Product Code:** RP-JR1-1

**Allergen:** rJug r 1 (*Juglans regia* allergen 1.0101)

**Lot No:** XXXXX

**Source:** *Pichia pastoris*

**Mol. Wt:** 15 kD

**Purification:** Purified from *Pichia pastoris* culture by multi-step chromatography. Purity > 95 % by SDS-PAGE.

**Concentration:** 1.2 mg/ml by OD 280; Extinction coefficient = 0.611.

**Formulation:** Preservative and carrier-free in 20 mM Tris, 0.3 M NaCl, pH 8.0. Filtered through 0.22µm filter.

**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** XXXXX

**Notes:** (1) rJug r 1 has a C-terminal 6xHis-tag.

**Gel Image:**

rJug r 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

[www.inbio.com](http://www.inbio.com)

700 Harris Street  
Charlottesville, VA  
22903  
United States

(434) 984-2304  
[mail@inbio.com](mailto:mail@inbio.com)

Vision Court  
Caxton Place  
Cardiff, Wales CF23 8HA  
United Kingdom

+44 (0) 2921 674 640  
[Info@inbio.com](mailto:Info@inbio.com)