

Product Code: NA-BG2-1

Allergen: Bla g 2 (Blattella germanica allergen 2)

Lot No: xxxxx

Source: B. germanica cockroach frass

Mol. Wt: ~36 kD

Purification: Purified from cockroach extract by monoclonal antibody affinity chromatography. Purity >95% by SDS PAGE.

Concentration: 2.0 mg/ml by OD 280; Extinction coefficient = 0.97.

Formulation: Preservative-free and carrier-free in phosphate buffered saline pH 7.4. Filtered through 0.22µ filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: xxxxx

Notes: (1) Natural Bla g 2 purifies as a 36 kD monomer and 72 kD dimer.

Gel Image:
nBla g 2



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com