

Product Code: NA-CF1-1

Allergen: nCan f 1 (*Canis familiaris* allergen 1)

Lot No: XXXXX

Source: Dog hair

Mol. Wt: 19 kD

Purification: Purified from dog extract by mAb affinity chromatography.
Purity > 95% by SDS-PAGE.

Concentration: 1.8 mg/ml by OD 280. Extinction coefficient = 0.75.

Formulation: Preservative and carrier-free in phosphate buffered saline, pH: 7.4.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: XXXXX

Notes: (1) Can f 1 belongs to the lipocalin protein family
(2) Not tested for endotoxin

Gel Image:

nCan f 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com