

Product Code: NA-DF1-1

Gel Image:

nDer f 1

Allergen: nDer f 1 (*Dermatophagoides farinae* allergen 1)

Lot No: XXXXX

Source: *D. farinae* culture medium

Mol. Wt: 24kD

Purification: From mite culture by affinity chromatography. Purity > 95 % by SDS-PAGE.

Concentration: 1.4 mg/ml by OD 280; Extinction coefficient = 1.9.

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µ filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: XXXXX

Notes: (1) Not tested for endotoxin



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com