

**Product Code:** NA-DF2-1

**Allergen:** Natural Der f 2  
(*Dermatophagoides farina* allergen 2)

**Lot No:** XXXXX

**Source:** *D. farinae* culture medium

**Mol. Wt:** 14kD

**Purification:** From mite culture by affinity chromatography.  
Purity > 95 % by SDS-PAGE.

**Concentration:** 1.1 mg/ml, by OD280 nm.  
Extinction coefficient = 0.6

**Formulation:** Preservative-free and carrier-free in phosphate buffered saline,  
pH 7.4. Filtered through 0.22µ filter.

**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** XXXXX

**Gel Image:**

nDer f 2



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

[www.inbio.com](http://www.inbio.com)

700 Harris Street  
Charlottesville, VA  
22903  
United States

(434) 984-2304  
[mail@inbio.com](mailto:mail@inbio.com)

Vision Court  
Caxton Place  
Cardiff, Wales CF23 8HA  
United Kingdom

+44 (0) 2921 674 640  
[Info@inbio.com](mailto:Info@inbio.com)