

Product Code: NA-DP2-1

Allergen: Natural Der p 2 (*Dermatophagoides pteronyssinus* allergen 2)

Lot No: XXXXX

Source: *D. pteronyssinus* culture medium

Mol. Wt: 14kD

Purification: From mite culture by affinity chromatography.
Purity > 95 % by SDS-PAGE.

Concentration: 1.8 mg/ml, based on OD 280 nm.
Extinction coefficient = 0.7

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µ filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: XXXXX

Gel Image:
nDer p 2



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com