

**Product Code:** NA-FD1-1

**Allergen:** nFel d 1 (*Felis domesticus* allergen 1)

**Lot No:** xxxxx

**Source:** cat hair

**Mol. Wt:** 18 kD

**Purification:** Purified from cat hair extract by affinity chromatography.  
Purity > 95% by SDS PAGE.

**Concentration:** 1.7 mg/ml by OD 280; Extinction coefficient = 0.435.

**Formulation:** Preservative-free and carrier-free in phosphate buffered saline,  
pH 7.4. Filtered through 0.22  $\mu$  filter.

**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** xxxxx

**Notes:** (1) Fel d 1 is glycosylated and migrates as a diffuse band  
on SDS-PAGE.

**Gel Image:**  
nFel d 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

[www.inbio.com](http://www.inbio.com)

700 Harris Street  
Charlottesville, VA  
22903  
United States

(434) 984-2304  
[mail@inbio.com](mailto:mail@inbio.com)

Vision Court  
Caxton Place  
Cardiff, Wales CF23 8HA  
United Kingdom

+44 (0) 2921 674 640  
[Info@inbio.com](mailto:Info@inbio.com)