

Product Code: RP-PA7-1

Gel Image:
rPer a 7

Allergen: rPer a 7.0102 (*Periplaneta americana* allergen 7)

Lot No: XXXXX

Source: *Pichia pastoris*

Mol. Wt: 33 kD; under reducing conditions

Purification: Purified from *Pichia pastoris* culture by multi-step chromatography.
Purity > 95% by SDS-PAGE.

Concentration: 0.73 mg/ml by BCA assay;
(Pierce™ BCA Protein Assay kit, Rockford, IL)

Formulation: Preservative-free and carrier-free in phosphate
buffered saline, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: XXXXX

Notes: (1) Made in USA.



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com