

Product Code: LTN-FD1-1

Allergen: nFel d 1 (*Felis domesticus* allergen 1)

Lot No: XXXXX

Source: cat hair

Mol. Wt: 18 kD

Purification: Purified from cat hair extract by affinity chromatography.
Purity > 95% by SDS PAGE.

Concentration: 1.3 mg/ml by OD 280; Extinction coefficient = 0.435.

Formulation: Preservative-free and carrier-free in phosphate buffered saline,
pH 7.4. Filtered through 0.22 µ filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: XXXXX

Notes: (1) Fel d 1 is glycosylated and migrates as a diffuse band
on SDS-PAGE.
(2) A LoTox™ product, Endotoxin < 0.03 EU/µg.

Gel Image:

LoTox™ nFel d 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com