

**Product Code:** NA-FD2-1

**Allergen:** nFel d 2 (*Felis domesticus* allergen 2)

**Lot No:** xxxxx

**Source:** Partially purified cat albumin

**Mol. Wt:** 69 kD

**Purification:** Purified by multi-step column chromatography.  
Purity >95% by SDS-PAGE.

**Concentration:** 1.1 mg/ml based on BCA Protein Assay.  
(Thermo Scientific, Rockford, IL)

**Formulation:** Preservative and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.


**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** xxxxx

**Notes:** (1) Fel d 2 is a serum albumin  
(2) Upper band is a Fel d 2 dimer  
(3) Not tested for endotoxin

**Gel Image:**

nFel d 2



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

[www.inbio.com](http://www.inbio.com)

700 Harris Street  
Charlottesville, VA  
22903  
United States

(434) 984-2304  
[mail@inbio.com](mailto:mail@inbio.com)

Vision Court  
Caxton Place  
Cardiff, Wales CF23 8HA  
United Kingdom

+44 (0) 2921 674 640  
[Info@inbio.com](mailto:Info@inbio.com)