

**Product Code:** LTN-DP2-1

**Allergen:** nDer p 2 (*Dermatophagoides pteronyssinus* allergen 2)

**Lot No:** XXXXX

**Source:** *D. pteronyssinus* culture medium

**Mol. Wt:** 14 kD

**Purification:** From mite culture by affinity chromatography.  
Purity > 95 % by SDS-PAGE.

**Concentration:** See product insert.

**Formulation:** Preservative-free and carrier-free in phosphate buffered saline,  
pH 7.4. Sterile filtered.

**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** See product insert.

**Notes:** (1) A LoTox™ product, Endotoxin < 0.03 EU/μg.

**Gel Image:**

LoTox™ nDer p 2



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

[www.inbio.com](http://www.inbio.com)

700 Harris Street  
Charlottesville, VA  
22903  
United States

(434) 984-2304  
[mail@inbio.com](mailto:mail@inbio.com)

Vision Court  
Caxton Place  
Cardiff, Wales CF23 8HA  
United Kingdom

+44 (0) 2921 674 640  
[Info@inbio.com](mailto:Info@inbio.com)