

Product Code: RP-DP2A-1

Allergen: rDer p 2 (*Dermatophagoides pteronyssinus* allergen 2.0101)

Lot No: XXXX

Source: *P. pastoris*

Mol. Wt: 14 kD

Purification: Purified from *Pichia pastoris* culture supernatant by multi-step chromatography. Purity > 95% by SDS-PAGE (See Panel).

Concentration: See product insert

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes:
(1) rDer p 2.0101 has a C-terminal 6xHis-tag.
(2) rDer p 2.0101 is suitable for T-cell or other in vitro studies, but does **not** react in the Group 2 mite ELISA.
(3) Not tested for endotoxin.

Gel Image:

rDer p 2



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com