

Product Code: AGAL-1

Allergen: Bovine Thyroglobulin (Galactose-alpha-1,3-galactose)

Lot No: XXXX

Source: Bovine thyroid gland.

Mol. Wt: ~300 kDa¹

Purification: Purity >95% as assessed by LC/MS-MS.

Concentration: See Product Insert


Formulation: Preservative and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See Product Insert

Notes: (1) Alpha-Gal forms tetramers with molecular weights of ~670 kDa.

Gel Image:



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com