

Product Code: BI-NAH2-1

Gel Image:

Biotin nAra h 2

Allergen: Biotinylated nAra h 2 (*Arachis hypogaea* allergen 2)

Lot No: XXXX

Source: Light roasted peanut flour (*Runner* cultivar)

Mol. Wt: 17-19 kD doublet

Purification: From delipidated peanut extract by multi-step chromatography.
Purity on SDS-PAGE >95%.

Concentration: See product insert.

Labelling Ratio: See product insert.

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4.
Sterile filtered.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert.

Notes: (1) Ara h 2 appears as a doublet ⁽⁴⁾.
(2) Avoid exposure to light.



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com