

**Product Code:** NA-AH2-1

**Allergen:** nAra h 2 (*Arachis hypogaea* allergen 2)

**Lot No:** XXXX

**Source:** Light roasted peanut flour (*Runner* cultivar)

**Mol. Wt:** 17-19 kD doublet

**Purification:** From delipidated peanut extract by multi-step chromatography. Purity >95% by SDS-PAGE.

**Concentration:** See product insert

**Formulation:** Preservative-free and carrier-free in phosphate buffered saline, pH 7.4. Sterile filtered.

**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** See product insert

**Notes:** (1) Ara h 2 appears as a doublet <sup>(4)</sup>.

**Gel Image:**

nAra h 2



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

[www.inbio.com](http://www.inbio.com)

700 Harris Street  
Charlottesville, VA  
22903  
United States

(434) 984-2304  
[mail@inbio.com](mailto:mail@inbio.com)

Vision Court  
Caxton Place  
Cardiff, Wales CF23 8HA  
United Kingdom

+44 (0) 2921 674 640  
[Info@inbio.com](mailto:Info@inbio.com)