

Product Code: NA-BE2-1

Allergen: nBer e 2 (*Bertholletia excelsa* allergen 2, Brazil nut)

Lot No: XXXX

Source: 11S globulin from Brazil nut.

Mol. Wt: Multiple bands ranging from 12 to 50 kD.

Purification: HPLC purified. Purity > 95 % by SDS-PAGE.

Concentration: See product insert.


Formulation: Preservative-free and carrier-free in, phosphate buffered saline, 1M NaCl, pH 7.4. Sterile filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert.

Notes: 11S globulins are synthesized as single peptide chain precursors that are cleaved post-translationally. Solubility of Ber e 2 is dependent on the presence of ≥1M NaCl.

Gel Image:
nBer e 2



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com