

Product Code: NA-BV1-1

Allergen: nBet v 1 (*Betula pendula* allergen 1)

Lot No: XXXX

Source: Birch pollen (*Betula pendula*)

Mol. Wt: 18 kDa

Purification: From birch pollen extract by multi-step chromatography. Purity > 95 % by SDS-PAGE

Concentration: 1.3 mg/ml by OD 280, extinction coefficient = 0.59.

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Notes: (1) Purified natural Bet v 1 appears as a triplet of 18-20kD proteins representing distinct isoforms of Bet v 1.

Gel Image:

nBet v 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com