

Product Code: NA-GAM1-1

Allergen: nGad m 1 (*Gadus morhua* allergen 1)

Lot No: XXXXX

Source: Atlantic cod

Mol. Wt: 11.45 kD

Purification: From Atlantic cod extract by salt precipitation and chromatography.
Purity > 95 % by SDS-PAGE.

Concentration: See product insert.

Formulation: Preservative-free and carrier-free in phosphate buffered
saline, pH 7.4. Sterile filtered.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Gel Image:

nGad m 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com