

Product Code: NA-PD6-1

Allergen: nPru du 6 (*Prunus dulcis* allergen 6)

Lot No: XXXX

Source: Almond flour

Mol. Wt: Multiple subunits ranging from 57-63kD.

Purification: Purified from almond flour by multi-step HPLC.
Purity > 95% by SDS-PAGE.


Concentration: See Product Insert

Formulation: Preservative and carrier-free in 20 mM Tris, 0.3M NaCl, pH 8.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: XXXX, when stored at -20°C.

Gel Image:



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com