

**Product Code:** RE-AG1-1

**Allergen:** rAln g 1 ( *Alnus glutinosa* allergen 1.0101, Alder)

**Lot No:** XXXX

**Source:** *E. coli*

**Mol. Wt:** 18 kD

**Purification:** Multi-step chromatographic purification. Purity >95% by SDS-PAGE.

**Concentration:** See product insert

**Formulation:** Preservative-free and carrier-free in 20mM Tris, 0.3M NaCl, 0.2M L-Arginine, pH 8.0. Filtered through 0.22µm filter.

**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** See product insert

**Notes:** (1) rAln g 1 has a C-terminal 6x His-tag.  
(2) rAln g 1 presents as a dimer under reducing conditions.

**Gel Image:**

rAln g 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

[www.inbio.com](http://www.inbio.com)

700 Harris Street  
Charlottesville, VA  
22903  
United States

(434) 984-2304  
[mail@inbio.com](mailto:mail@inbio.com)

Vision Court  
Caxton Place  
Cardiff, Wales CF23 8HA  
United Kingdom

+44 (0) 2921 674 640  
[Info@inbio.com](mailto:Info@inbio.com)