

Product Code: RE-CF4-1

Allergen: rCan f 4 (*Canis familiaris* allergen 4.0101)

Lot No: xxxx

Source: *E. coli*

Mol. Wt: 18.5 kD

Purification: Multi-step chromatographic purification.
Purity > 95% by SDS-PAGE.

Concentration: See product insert

Formulation: Preservative and carrier-free in phosphate buffered saline, pH 7.4.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: 1. rCan f 4 has a N-Terminus 6XHis-tag.

Gel Image:

rCan f 4



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com