

Product Code: RE-JR5-1

Allergen: rJug r 5 (*Juglans regia* Allergen 5.0101)

Lot No: XXXX

Source: *E. coli*

Mol. Wt: 19 kD

Gel Image:
rJug r 5



Purification: Multi-step chromatographic purification.
Purity > 95 % by SDS-PAGE.

Concentration: See Product Insert

Formulation: Preservative-free and carrier-free in 0.2M L-Arginine, phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See Product Insert

Notes: (1) rJug r 5 has a C-terminal 6xHis-tag.

For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com