

Product Code: RP-AG1-1

Allergen: rApi g 1 (*Apium graveolens*, allergen 1.0101)

Lot No: XXXX

Source: *P. pastoris*

Mol. Wt: ~16 kD

Purification: Purified from *Pichia pastoris* culture by multi-step chromatography. Purity > 95 % by SDS-PAGE, run under reducing conditions.

Concentration: See product insert

Formulation: Preservative and carrier free in 1X PBS buffer, pH 7.4
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: (1) rApi g 1 has a C-terminal 6xHis-tag.

Gel Image:

rApi g 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com