

Product Code: RP-DF2-1

Allergen: rDer f 2 (*Dermatophagoides farinae* allergen 2)

Lot No: XXXX

Source: *P. pastoris*

Mol. Wt: 14 kD

Purification: Purified from *Pichia pastoris* culture supernatant by affinity chromatography. Purity > 95% by SDS-PAGE.

Concentration: See product insert

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: (1) Not tested for endotoxin.

Gel Image:

rDer f 2



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com