

Product Code: RP-DP7-1

Allergen: rDer p 7 (*Dermatophagoides pteronyssinus* allergen 7)

Lot No: XXXX

Source: *P. pastoris*

Mol. Wt: 22kD and 27kD

Purification: From *Pichia pastoris* culture by two step column chromatography.
Purity > 95% by SDS-PAGE (See Panel).

Concentration: See product insert

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: (1) rDer p 7 is produced in *P. pastoris* as two bands: at 22kD and 27kD.

Gel Image:

rDer p 7



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com