

Product Code: RP-EM1-1

Allergen: rEur m 1 (*Euroglyphus maynei* allergen 1.0101)

Lot No: XXXX

Source: *P. pastoris*
(Clone N16Q, N133Q, N-glycosylation sites mutated)

Mol. Wt: 36 kD

Purification: Purified from *Pichia pastoris* culture supernatant by multi- step chromatography. Purity > 95% by SDS-PAGE (See Panel).

Concentration: See Product Insert

Formulation: Preservative-free and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See Product Insert

Notes: (1) rEur m 1 has a C-terminal 6xHis-tag.
(2) N-glycosylation motifs at positions 16 and 133 were mutated.

Gel Image:
rEur m 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com