

**Product Code:** RP-GM4-1

**Allergen:** rGly m 4 (*Glycine max* allergen 4.0101)

**Lot No:** XXXX

**Source:** Pathogenesis-related protein expressed in *Pichia pastoris*

**Mol. Wt:** 17 kD

**Purification:** Purified from *Pichia pastoris* by multi-step HPLC. Purity > 95 % by SDS-PAGE.

**Concentration:** See product insert

**Formulation:** Preservative and carrier-free in PBS, pH 7.4. Filtered through 0.22µm filter.

**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** See product insert

**Notes:** (1) rGly m 4 has a C-terminal 8xHis-tag.

**Gel Image:**

rGly m 4



**For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.**

An InBio® product. Made in the USA

[www.inbio.com](http://www.inbio.com)

700 Harris Street  
Charlottesville, VA  
22903  
United States

(434) 984-2304  
[mail@inbio.com](mailto:mail@inbio.com)

Vision Court  
Caxton Place  
Cardiff, Wales CF23 8HA  
United Kingdom

+44 (0) 2921 674 640  
[Info@inbio.com](mailto:Info@inbio.com)