

## Recombinant Jug r 1

**Product Code:** RP-JR1-1

**Allergen:** rJug r 1 (*Juglans regia* allergen 1.0101)

**Lot No:** XXXX

**Source:** *Pichia pastoris*

**Mol. Wt:** 15 kD


**Purification:** Purified from *Pichia pastoris* culture by multi-step chromatography. Purity > 95 % by SDS-PAGE.

**Concentration:** See product insert

**Formulation:** Preservative and carrier-free in 20 mM Tris, 0.3 M NaCl, pH 8.0. Filtered through 0.22µm filter.

**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** See product insert

**Released by:** 

**Notes:** rJug r 1 has a C-terminal 6xHis-tag.

**Gel Image:**  
rJug r 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

[www.inbio.com](http://www.inbio.com)

700 Harris Street  
Charlottesville, VA  
22903  
United States  
(434) 984-2304  
[mail@inbio.com](mailto:mail@inbio.com)

Vision Court  
Caxton Place  
Cardiff, Wales CF23 8HA  
United Kingdom  
+44 (0) 2921 674 640  
[Info@inbio.com](mailto:Info@inbio.com)