

Product Code: RPI-PP1-1

Allergen: rPin p 1 (*Pinus pinea* allergen 1.0101)

Lot No: XXXX

Source: *Pichia pastoris*

Mol. Wt: 17.4 kD

Purification: Purified from *Pichia* culture supernatant using multi-step chromatography. Purity >95% by SDS-PAGE.

Concentration: See product insert

Formulation: 20mM Tris-HCl, pH 8.5. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: (1) rPin p 1 has a N-terminal 6xHis-tag.

Gel Image:

rPin p 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com