

## LoTox™ Egg-White Protein

**Product Code:** LTN-GDE-1


**Allergen:** *Gallus domesticus (G. gallus)*

**Lot No:** XXXXX

**Source:** Natural egg-white powder

**Mol. Wt:** Gal d 1: 28kD, Gal d 2: 44kD, Gal d 3: 78kD, Gal d 4: 14kD

### Gel Image



Egg-white Protein

**LC-MS/MS:** See product insert.

**ELISA:** See product insert.

**Concentration:** See product insert.

**Formulation:** Preservative and carrier-free in sterile, 0.5M NaCl, 0.1M Potassium Phosphate pH 7.40. Sterile filtered.

**Storage:** Store at -20°C. Re-freeze unused portion and aliquot into smaller quantities to avoid freeze/thaw cycles. Gently defrost at room temperature.

**Expiration:** See product insert.

**Released:**

*Liz D Vaile*

**Notes:**  
(1) (1) A LoTox™ product: Endotoxin  $\leq 0.03$  EU/ $\mu$ g.  
(2) Relative abundance results are based on label-free quantification of precursor peptide ion intensity (peak area, incl. unique + razor peptides) using Proteome Discoverer 2.2 algorithms.  
(3) Calculations of relative abundance can be variable and depend on multiple factors such as sample complexity, instrument resolution and analysis software.

For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

[www.inbio.com](http://www.inbio.com)

700 Harris Street  
Charlottesville, VA 22903  
United States  
(434) 984-2304  
[mail@inbio.com](mailto:mail@inbio.com)

Vision Court, Caxton Place  
Cardiff, Wales CF23 8HA  
United Kingdom  
+44 (0) 2921 674 640  
[Info@inbio.com](mailto:Info@inbio.com)