

**Product Code:** RPI-FA1-1

**Allergen:** rFra a 1 (*Fragaria ananassa* allergen 1.0102)

**Lot No:** XXXX

**Source:** *Pichia pastoris*

**Mol. Wt:** 19 kD

**Purification:** Purified from *Pichia* culture supernatant using multi-step chromatography. Purity >95% by SDS-PAGE.

**Concentration:** mg/ml by. Extinction coefficient = 0.773

**Formulation:** 20mM Tris, 0.3M NaCl, pH 8.0. Filtered through 0.22µm filter.

**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** ; when stored at -20°C.

**Notes:** (1) rFra a 1 has a N-terminal 6xHis-tag.

**Gel Image:**

rFra a 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

[www.inbio.com](http://www.inbio.com)

700 Harris Street  
Charlottesville, VA  
22903  
United States

(434) 984-2304  
[mail@inbio.com](mailto:mail@inbio.com)

Vision Court  
Caxton Place  
Cardiff, Wales CF23 8HA  
United Kingdom

+44 (0) 2921 674 640  
[Info@inbio.com](mailto:Info@inbio.com)