

Product Code: RE-AM1-1

Allergen: rApi m 1 (*Apis mellifera* allergen 1.0101)

Lot No: 43243

Source: *E. coli*

Mol. Wt: 16 kD

Purification: Purified from *E. coli* culture by multi-step HPLC chromatography. Purity > 95 % by SDS-PAGE.

Concentration: 0.3 mg/ml by OD 280; Extinction coefficient = 1.45.

Formulation: Preservative and carrier-free in 20mM Tris, 0.3M NaCl, 0.2M L-arginine pH 8.0. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: 03/2029; when stored at -20°C.

Notes: (1) rApi m 1 has a C-terminal 6xHis-tag.

Gel Image:

rApi m 1



For research and commercial use *in vitro*: not for human *in vivo* or therapeutic use.

An InBio® product. Made in the USA

www.inbio.com

700 Harris Street
Charlottesville, VA
22903
United States

(434) 984-2304
mail@inbio.com

Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

+44 (0) 2921 674 640
Info@inbio.com